

### ▶ Pradhan Karma

We have quality products in Pradhan Karma. All our models are efficient and display quality in Panchakarma (Ayurvedic) products.



#### **IMI - 2256: VAMAN PEETH** (Chair and Table)



Tubular steel Vaman table is mounted on 4 wheels & fitted with stainless steel sink, side trays and comes with measuring Jar for vaman kriya.

A steel chair with wheels and reclining back is also provided for comfort.

[\[Request Quotation / More info...\]](#)

#### **IMI - 2258: NASYA PEETH**



Tubular steel chair with wooden rungs seat and back is comfortable for the patient to sit in semi-recumbent position and place his head on the adjustable head support for the easyness of the doctor to put medicine in the nostrils.

#### FEATURES:

\*Frame: Constructed of square tube and finished oven baked enamel.

\*Seat & Back: Well polished Teak-wood rungs fitted for comfort, back

has adjustable recline & also fitted with Head support.

[\[Request Quotation / More info...\]](#) 

---

### **IMI - 2259: VEERACHAN CHAIR**



Comfortable commode chair, constructed of square pipe and finished in oven baked enamel, is fitted with solid Teak wood seat with plastic commode and back with adjustable recline, is fitted with teak wood rungs, the side-arm supports adds comfort to the patient.

Chair is mounted on 10cm dia. heavy duty casters for mobility.

[\[Request Quotation / More info...\]](#) 