

Rehabilitation Aids

▶ Walking Sticks / Canes

Presenting walking sticks / canes in Aluminium with anodized bright polish. All models are strong, sturdy and fitted with 'L' shaped handle for better grip. India Medico's walking sticks / canes have good demand in India and also exported to countries all over the world.



IMI - 3061: Walking Cane / Walking Stick (Aluminium, Fixed)



Anodized aluminium bright polished fixed height cane. Fitted with L shape handle for better grip.

[\[Request Quotation / More info...\]](#) 🏠

IMI - 3062: Walking Cane / Walking Stick (Aluminium, Adjustable)



Adjustable in height by slipping spring locking pin. Fitted with PVC tip. Anodized bright polish finish, with wooden handle for better grip (also available in 'U' shape and 'L' shape handle).

Packing Details:

CARTON: 36" x 24" x 22" for 80 pieces.

[\[Request Quotation / More info...\]](#) 📄

IMI - 3065: Walking Cane / Walking Stick (Four legged)



Metallic stick with "U" shape hand grip. Broad base with specially angled legs, can be held close to user body for maximum support. Height is adjustable from 75 c.m. to 95 c.m.

[\[Request Quotation / More info...\]](#) 

IMI - 3070: Tips for Walking Cane / Crutch (19 mm, 22 mm, 25 mm diameter)



Broad based non- slip moulded tip, reinforced with a steel insert that is moulded into the base. Suitable for cane and crutches. Available in 3 internal sizes of 19 mm, 22 mm & 25 mm.

[\[Request Quotation / More info...\]](#) 